



Salient Features

- Auto Calibration with External Weight
- Easy to read LCD Display with backlight
- Multifunction Weighing Units
- Auto Zero tracking, Low Battery Indicator
- Overload Protection Design
- S.S. Weighing Pan, Level indicator, Adjustable Fit
- High Resolution, Quick Weighing, Accurate Result
- In built rechargeable battery for continues use.

Precision Balance

PGB3110 & 6100

HSN Code - 9016

HSN Code - 8423

SPECIFICATIONS	MODEL	PGB3110	PGB6100
	Capacity	3000 g	6000 g
	Readability	0.01 g (10 mg)	0.1 g (100 mg)
	Repeatability (±)	0.01 g	0.1 g
	Linearity (±)	0.02 g	0.2 g
	Pan Size	160 x 160 mm	
	Make	WENSAR	
	Display	LCD	
	Calibration	Automatic External	
	Tare Range	Full	
	Operating Temperature	10°C to 40°C	
	Power supply	AC 220V/ 50-60 Hz	
	Dimension (LxWxH) mm	440 x 265 x 220	
	Weight	4 kg	

★ Specifications are subject to change without notice.